

Identification and quantification of preservative chemicals in common household products

Session 3

The aim of this session is to analyse the extract that you purified during Session 2. You will identify and quantify the parabens present in your extract and therefore be in a position to assess whether the product you have decided to investigate complies with the 76/768/EEC council directive relating to the use of preservatives in cosmetic products.

Using your findings from the Session 3 pre-laboratory exercise, discuss the aims of the current session and suggest some methods that should enable you to achieve these aims. You will have **15 minutes** to do this.

Note: When designing your experimental procedure, consider the following:

- The desired outcome(s)
- Glassware/equipment/solvents that are available to you
- Time
- Safety
- Data/information to be recorded

Identification and quantification of preservative chemicals in common household products

Method to carry out analysis of extracts containing parabens

Parabens can be analysed using GC or GC-MS techniques. However, they need to be derivatised prior to the analysis. In the case of parabens, a suitable silylation (derivatising) agent is bis-trimethylsilyl trifluoroacetamide (MTBSTFA).

Silylation procedure:

1. Add 10 µg of each of the paraben standards to a 2 mL GC vial (use the solution that you prepared during Session 2). Evaporate the solvent and add 1 mL of ethyl acetate.
2. Dissolve your extract in 1 mL of ethyl acetate and transfer 10 µL of the solution into a 2 mL GC vial. Evaporate the solvent and add exactly 1 mL of ethyl acetate to the vial.
3. Add 1 mL of ethyl acetate to your procedural blank extract.
4. Add 40 µL of MTBSTFA into each vial and heat them at 70°C for 30 minutes in a heating block in a fume cupboard. When cool, your 3 samples will be ready for analysis by GC or GC-MS.

Identification and quantification of preservative chemicals in common household products

Author	Simon Belt, Guillaume Massé
Title	Parabens in household products
Classification	Laboratory Manuals - Chemistry
Keywords	ukoer, parabens, chemistry, analytical, GC
Description	Individual lab sheets - Tutor
Creative Commons Licence (url)	http://creativecommons.org/licenses/by-nc-sa/2.0/uk/
Language	English
File size	120 kB
File format	pdf